

Thompson Retractor

DO'S AND DON'TS OF INSTRUMENT CLEANING AND STORAGE

DO

- ✓ Check that all parts have been properly cleaned and sterilized and that no unacceptable stains or residues are present on parts.
- ✓ Always pre-clean instruments before placing into sterilizer.
- ✓ Use low to average pH products/cleaning agents.
- ✓ Use distilled water when soaking instruments.
- ✓ Always ensure instruments are completely dry before storing.
- ✓ Always store instruments in clean condition; in a dry & sterile environment.
- ✓ Protect instruments when storing; wrap instruments or use instrument case trays.
- ✓ Always store instrument handles UNLOCKED (open) when not in use.
- ✓ Always sterilize instrument handles with cam joint UNLOCKED (open).
- ✓ Always UNLOCK (open) handles before adjusting rods or disassembling equipment.
- ✓ Ensure cam joint holes have smooth inner surface, free of scratches or burs, before inserting rods.
- ✓ Ensure rod surfaces are smooth, free of gouges/scoring, before placing into cam joint holes.
- ✓ Carefully align screw knobs to smoothly tighten/engage threads. (Keep threads clean and lubricated.)
- ✓ Always depress blade handle release button while attaching/detaching blades.

DO NOT

- ✗ Use parts with residue—call Thompson if you have questions on staining or residues.
- ✗ Sterilize instruments without pre-cleaning; blood/debris can embed and cross-contaminate.
- ✗ Use high pH products or those containing chloride; will strip protective layer on metal and cause rust.
- ✗ Use saline when soaking instruments.
- ✗ Store instruments wet or semi-dry; will stain and corrode the metal.
- ✗ Store instruments dirty; debris can contaminate and rust the metal.
- ✗ Store instrument handles LOCKED (closed) when not in use; will damage/collapse cam joints.
- ✗ Sterilize instrument handles with cam joint LOCKED (closed); will damage/collapse cam joints.
- ✗ Adjust rods while cam joint handles are LOCKED down; twisting/force will damage both pieces.
- ✗ Insert rods into damaged cam joint holes (force fit); will damage both pieces.
- ✗ Insert damaged rods into cam joint holes; will damage both pieces.
- ✗ Screw down knobs with threads misaligned; will strip/cross threads (dirty threads will jam knobs).
- ✗ Pull blades from handles without depressing handle release button; will damage blade & nipple.
- ✗ Use tools (i.e. screwdriver) to disassemble instruments for cleaning.

 **Thompson**
SURGICAL INSTRUMENTS
10341 East Cherry Bend Road
Traverse City, Michigan 49684
phone: 231.922.0177
fax: 231.922.0174
thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2020 Thompson Surgical Instruments, Inc.

® S-Lock®, ® PLA®, and the Thompson "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.
Patents: trpat.com

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
THE NETHERLANDS

CE0297

Rev F
071620
dd

**Made In
USA**