आगे का कवर

अनवािरय अनावरण

Thompson तकनीक: ब्रु (रविर्स लपि) अग्रमि कमर पहुंच

फ़ोटो इनसेट: कशेरुक शरीर के पार्श्व पहलू को संलग्न कर के लिप स्थरिता को खीचे।

"Thompson AL Brau ब्लेड सुरक्षा के 15 वर्ष ट्रैक रिकॉर्ड के साथ कमरदर्द रीढ़ के अग्रिम छोटे-खुले दुष्टिकोण के लिए आवश्यक हो गए है। आपको स्थिर, ठोस शीला जैसा अनावरण प्राप्त होता है जिससे 'फडिल फैक्टर' को कम करके इस तरह से समय को फिर से समायोजित करने की आवश्यकता के बिना उत्कृष्ट दुश्य उपलब्ध कराया जाता है। "

- सल्वाडोर ए. ब्रू, एम.डी, एफ.ए.सी.एस.

#### पषठ 2:

THOMPSON तकनीक

सल्वाडोर ए. ब्रु., एम.डी, एफ.ए.सी.एस.

हों. बरू न्यूयॉर्क शहर के माउंट सिनाई अस्पताल में सामान्य संवहनी की सर्जरी में प्रशक्षित है और 20 से अधिक वर्षों तक रीढ़ की हहड़ी के सभी सुतरों के लिए अग्रिम पहुंच प्रक्रियों को प्रदर्शति कर रहे हैं। उन्होंने कमर की रीढ़ की हहड़ी के लिए 1500 से अधिक खुले अन्तर्नाष्ठ दुष्टिकोण का प्रदर्शन किया है और अग्रिम लैप्रोस्कोपिक और धोरैकोस्कोपिक पहुंच में भी अग्रणी रहे हैं। अब सेवानविृत्त, डॉ. ब्रू यू.एस.सी. में सर्जरी के सहायक नैदानिक प्रोफेसर थे, और लॉस एंजल्सि के सिंडर-सिनाई मेडिकल सेटर में काम किया है।

डॉ. ब्रू ने कमर की रीढ़ की हहडी के लिए एक नया "छोटा-खुला" अग्रिम दुष्टिकोण विकसित किया है। 2000 के अक्टूबर में उन्होंने न्यू ऑरलियन्स में नॉर्थ अमेरिकन स्पिन सोसाइटी की वार्षिक बैठक में 386 मामलों में इस दुष्टिकोण के साथ अपने अनुभव को वीडियी के साथ प्रस्तुत किया। यह प्रक्रिया डॉ. ब्रू द्वारा रूपांकित किए गए विशेष पुनरावृत्तियी पर निर्भर है। इस प्रत्याकर्षक ने 9 अक्टूबर 2001 को पेटेट प्राप्त किया, और Thompson Surgical Instruments, Inc. से विशेष रूप से उपलब्ध है।

कमर की रीढ़ की हड्डी का दुष्टिकोण मुख्य सर्जन की अनावरण जल्दी से प्रदान करने की क्षमता पर और सुरक्षित रूप से संवहनी क्षति की घटना की रिपीर्ट पर निर्भर करता है। पिजरों और कृत्रिम इसिक्स के सेरेखण के लिए "सीधा" आगे-पीछे अनावरण की आवश्यकता ने एक छोटा चीरा प्रदान करने के लिए मुख्य सर्जन के लिए एक महत्वपूरण चुनौती पेश की है और अभी तक श्रोणिफलकीय पात्र और स्वायत्त तंत्रिका जाल की चोट को रोकने के लिए आवश्यक सुरक्षा की डिग्रि को बनाए रखना है। Thompson प्रत्याकर्षक प्रणाली के उपयोग में वर्णित दृष्टिकोण ने इन चिताओं को काफी कम कर दिया है।

वशिष रविरस लपि बुलेड

कशेरुक शरीर के पार्श्व पहलु को संलग्न कर के लिप सुथरिता को खीचे।

अनवारय अनावरण

रेडयिलिसिट रविर्स लिप अग्रिम कमर ब्लेड पूर्वकाल के कमर की रीढ़ की अभूतपूर्व अनावरण की पेशकश करता है और ब्लेड की फसिलन को रोकता है।

रेडियोल्यूसेन्ट, सशक्त ब्लेड

बेहतर शक्ति भारी खिचाव के तहत अत्यथिक झुकाव को रोकती है। सभी रोगी आकारों को समायोजित करने के लिए 250 मिमी तक की गहराई तक उपलबप है।

टप्पणी:आदेश की जानकारी देने के लिए पृष्ठ 7 देखें।

टपि्पणी:क्योकि हम लगातार सर्वोत्तम उत्पादों को उपलब्ध कराने का प्रयास करते हैं, इस उपयोगकर्ता पुस्तिका में द्विए कुछ चित्र प्राप्त उत्पाद से थोड़े भन्नि दक्षिाई दे सकते हैं।

THOMPSON तकनीक:

प्रश्न या अतरिक्त जानकारी:1.800.227.7543

पषठ 7:

रटि्रोपीरटिोनयिल छोटे-दुरष्टिकोण सेटअप चरण

नीचे डॉ. ब्रू द्वारा उल्लखिति, रट्रिरोपीरटिोनयिल अग्रमि कमर के अनावरण के लिए सुझाए गए कदम उठाए गए है।

चरण 1

रोगी को एक एक्सरे टेबल पर चित्त स्थिति में रखें।

दुष्टिकीण सर्जन बाई तरफ खड़ा है और दाई ओर सहायक है। क्रनिओकौडाड पुलेन में अनुप्रस्थ चीरा का सुतर संपर्क होने वाली रीढ़ की हड्डी के सुतर पर निर्भर करता है। इस चीरा का उचति सुधान निर्धारति करने के लिए रीढ़ की एक पार्श्व एक्सरे आवश्यक है।

सावधानी:यद्िरोगी मोटा है, तो रेल क्लंप रखने पर अहरि तंत्रिका को दुबाने से बचें।

टिप: जब आवश्यक हो, हमारे रेल एक्सटेशन (# 41 9 17) का उपयोग करके वाइडर या तालिका का उपयोग करे या तालिका की चौड़ाई में 2 ¾ "जोडे।

चरण 2

बाई रेक्टस मांसपेशायी को गतिशील रूप से जुटाया जाता है। रेक्टस मांसपेशायी को शुरू में मध्यवर्ती में वापस ले कर, जब तक का पैरिटीनियम को चमकते हुए देखा जाए, तब तक 4 से 5 सेटीमीटर ट्रांस्वार्सलिस प्रावरणी के पीछे वाली कंधे को ध्यान में रखते है। एक हेमोस्टेट के साथ किनारों को पकड़े और इसे उठाएं और पेरिटीनियम से हो तो बहुत सावधानी से विच्छेदन करें। जितना संभव हो उतना कमजोर और श्रेष्ठ रूप से इसे अलग कर लें। अपनी तर्जनी का उपयोग कर के, सावधानीपूरवक फेसिअल चीरे के किनारे पर पेरिटीनियम को पीछे से धक्का दे इसके और आंतरिक तरिछा के निचले हिस्से और ट्रांसवरसस मांसपेशायी और प्रावरणी बीच एक सतह तैयार करें। यह आपको रेट्रोपरिटीनियल स्पेस पर ले जाएगा।

बाद में सावधानीपुर्वक कुंद उंगली विच्छेदन को जारी रखें, और फिर पेरितीनीयम को नितम्ब की मस्पेशियों से दूर हटाने बीच से धक्का देना शुरू करें। सावधान रहें कि इस समय पर रिट्रोपसोस क्षेत्र में प्रवेश न करें, क्योंकि इससे अंधा थैली में अनावश्यक रक्तस्त्राव हो सकता है। नितम्ब की मस्पेशियों के ऊपर जेनीटोफेमोरल तंत्रिका को आसानी से पहचाना जा सकता है। पेरितीनीयम को नितम्ब की मस्पेशियों से दूर हटाते समय मृत्रनली को पहचाना जा सकता है। इन दोनों संरचनाओं को क्षति से संस्कृषित किया जाना चाहिए।

फ़ोटो इनससेट: रेट्रोपेरटिोनीयल वचिछेदन का रेक्टस मांसपेशी के पार्श्व से प्रारंभ करें

thompsonsurgical.com पर हमारी ऑनलाइन मुलाकात लें

#### पषठ ४

THOMPSON तकनीक/प्रत्याकर्षक ब्लेड लगाना रेट्रोपेरिटीनीयल लघु-उपागम सेट अप के चरण (जारी)

#### वरण ३

नितम्ब की मस्पेशियों की पहचान होने पर, इिस्क और कशेरुकाय शरीर और एलिएक धमनी को महसूस करने के लिए मध्यवर्ती स्पर्श करे। इस समय पर, यहि चीरे का आकार अनुमति देता है, तो रेट्येपरिटीनियल क्षेत्र में पूरा हाथ डाले और मुद्री बनाएं। बंदिस्त मुद्री के ऊपर और नीचे से सभी दिशाओं में पैरिटीनियम दूर करने के लिए घुमाएं। आम और बाह्य एलिएक धमनियों की पूरी लंबाई को जितना हो सके उतनी दूर तक खोलने के लिए कुंद विच्छेदन जारी रखे, और फिर धमनी के पार्श्व किनारे पर सावधानीपूरवक कुंद विच्छेदन जारी रखे। सामान्य एलिएक नस को खोलेगा। एलीयोलम्बर शारी (री) की पहचान करने के लिए बाद में विच्छेदन जारी रखे। सामान्य एलिएक नस और लम्बर की नसी के गठन में भिन्ताएं आम होती है, और ये शारीओं की पहचान, लैगेट करने और आड़ी काटने करने और अवसाद से बचने के लिए बहुत सावधानी बरतनी चाहिए। बाएं एलिएक नस और धमनी को पीनट सूर्पज, उंगलियों और कुंद लिएट विच्छेदन का उपयोग करके रीढ़ से अलग किया जा सकता है।

फोटो इन्सेट: रेट्रोपेरीटोनयिल दुष्टिकोण जारी

#### नरमा ४

जिस तरफ़ शल्प-चिकित्सक का परिचालन क्षेत्र खाली रखता है उस तरफ जीवाणुरहृति कपड़े पर टेबल रेल पर एलटि रेल क्लैप रखें।(ए)। क्रॉसबार को ऑपरेटिव साइट के 2 cm ऊपर जॉइंट में डाला जाता है(बी)।

#### चरण व

सभी संवहनी संरचनाओं को फरि बाएं से दाहिनी ओर हिलाया जाता है, जो शामिल डिस्क और करोरुका का पर्याप्त दूश्य प्रदान करता है। निकायों के पूर्वकाल सतह पर वेळी में पुमने वाले हिस्सों वाले वेसल क्लिप के बीच की तरफ जा सकता है और कुंद विच्छेदन के साथ पक्षों में फिसल सकता है। सुनिश्चित करें कि आप शरिं। और लिग्जिमेंट के बीच कम से कम एक उंगली पा सकते हैं ताकि आप सही पार्श्व को छू सके। आपकी उंगली से ऊपर वाले वेसल के साथ रीढ़ की हड़डी के किनारे रख सकते हैं।

फोटो इन्सेट: एलीयो लम्बर नस के बंधाव से बाद ट्यूब के नीचे ऊँगली का विच्छेदन

THOMPSON तकनीक: अग्रीम लम्बर /सल्वाडोर ए. ब्र्, एम.डी, एफ.ए.सी.एस. प्रश्न या अतरिकृत जानकारी:1.800.227.7543

#### पषठ 5

रटिरोपीरटीनियल छोटे-दुरष्टिकोण सेटअप चरण (जारी)

#### वरण 6

पार्श्व वस्तिर हाथ क्रॉसबार (ए) से जोड़े जाते है और रोगी के क्षतिजि (बी) के ऊपर रखें जाते है।

#### चरण 7

रेक्टस को पाशर्ववत करने के बाद सर्जन का बाया हाथ फरि से रेट्रोपेरिटीनयिल क्षेत्र में फरि से प्रवेशकरेगा, और उंगलियां रीढ़ की हहड़ी के दाहानी ओर अपना रास्ता ढूंढेगी। एक रेडियीलिसेट रविर्स लिप अग्रिम रीढ़ के ऐक्सेस ब्लेड को रीढ़ की हहड़ी के दाई तरफ एक उंगली का मारगदुरशक के रूप में उपयोग करके रखा जाएगा।

फोटो इनसेट: रीढ़ की हड्डी के दाहिनी ओर प्रारंभिक रविर्स लिप ब्लेड की प्रविष्ट

### चरण 8

इस बुलेंड को फारे प्रत्याकर्षक फ्रेम के पार्शव विस्तार हाथों से जोड़ा जाता है, संवहनी सरचना को ऊपर उठाना और रीढ़ की हड़ड़ी की सतह को खोलता है। Thompson प्रत्याकर्षक से सुरक्षित तरीके से जोड़ने के बाद, रविर्स-लिप किए गए बुलेंड हिलेंगे नहीं। रविरस लिप बुलेंड के किनोरे को रीढ़ की हड़ड़ी के किनोरे पर लगा के रखता है और तनाव देने पर फिसलने से रोकता है। इस रविर्स लिप के बिना, प्रत्याकर्षक ब्लेंड प्रभावी ढंग से काम नहीं करेगा।

### चरण 9

रेक्ट्स को पाशर्ववत पीछे हटने के बाद, सीधा ए.पी दुश्य में रीढ़ की हहड़ी को खोलने के लिए प्रत्याकर्षक ब्लेड को धक्का देते समय कम प्रतिरिध होगा और एक धूरेडे किए हुए उपकरण, फेमोरल रिग या कृत्रिम इसिक को सम्मिलित करने के लिए स्लीव को रखने की अनुमत दिते हैं। रीढ़ की बाई ओर एक दूसरा रिवर्स लिप ब्लेड रखे और Thompson फ्रेंम से जोड़े। सामान्यतया, अतिरिक्त प्रत्याकर्षक ब्लेड को अनावरण पूरा करने के लिए श्रेष्ठ और / या न्यून रूप से रखा जाना चाहिए। With the blades well anchored to the lateral wall of the vertebral column, the spine surgeon and the assistant can now work on the disc without other hands or retractors being in the way and with relative security that vessels will not move around the retractors and expose themselves to injury.

फ़ोटो इनसेट: दोनो रविर्स लिप ब्लेड तैनात किए गए और डिस्क के पार्श्व पहलू और कशेरुक शरीर पर लगे हुए है।

thompsonsurgical.com पर हमारी ऑनलाइन मुलाकात लें

पृष्ठ 6 + 7

THOMPSON तकनीक / अग्रमि कमर प्रणाली + वैकल्पिक सहायक उपकरण आदेश की जानकारी


पीछे का कवर

परतीक कविदंती:

प्रतीक कविदंती: नर्माता अधिकृत EC Rep CE Mark Mark चेतावनयिँ / सावधानयिँ गैर-जीवाणुरहति



**Uncompromised Exposure**<sup>™</sup>



# **Thompson Techniques:**

# **BRAU (REVERSE LIP) ANTERIOR LUMBAR ACCESS**

"The Thompson AL Brau Blades have become essential for the Anterior Mini-Open approach to the lumbar spine with a proven 15 year track record of safety. You get steady, rock solid exposure providing excellent visualization throughout without the need to re-adjust thus saving time by reducing the 'fiddle factor.'"




# Salvador A. Brau, MD, FACS

Dr. Brau trained in general vascular surgery at the Mount Sinai Hospital in New York City and has been performing anterior access procedures to all levels of the spine for over 20 years. He has performed over 1,500 open anterior approaches to the lumbar spine and has also been a pioneer in anterior laparoscopic and thoracoscopic access. Now retired, Dr. Brau was Assistant Clinical Professor of Surgery at USC, and worked at Cedars-Sinai Medical Center in Los Angeles.

Dr. Brau has developed a new "mini-open" anterior approach to the lumbar spine. In October of 2000, he presented his experience with this approach in 386 cases, along with a video of the procedure, to the North American Spine Society Annual Meeting in New Orleans. This procedure is dependent on special retractors designed by Dr. Brau. These retractors received a patent on October 9, 2001, and are available exclusively from Thompson Surgical Instruments, Inc.

The anterior approach to the lumbar spine is heavily dependent on the ability of the access surgeon to provide exposure quickly and safely in view of a reported incidence of vascular injury. The requirement of a "straight on" anterior-posterior exposure for alignment of cages and artificial discs has presented a significant challenge for the approach surgeon to provide a small incision and yet maintain the degree of safety necessary to prevent injury to the iliac vessels and autonomic nerve plexus. The approach described here utilizing the Thompson retractor system significantly reduces these concerns.



**Exclusive Reverse Lip Blades** 

Reverse lip stabilizes retraction by engaging the lateral aspect of the vertebral body.



**Uncompromised Exposure** 

Radiolucent Reverse Lip Anterior Lumbar blades offer phenomenal exposure of the anterior lumbar spine and prevent blade slippage.



### Radiolucent, Strong Blades

Superior strength prevents excessive flexing under heavy retraction.

Available up to 250mm deep to accommodate all patient sizes.



### NOTE

See page 7 for ordering information.



### NOTE

As we continually strive to provide the best products possible, some of the images in this user manual may appear slightly different from the product received.

# Retroperitoneal Mini-Approach Set Up Steps

Below are the suggested set up steps for retroperitoneal anterior lumbar exposure, as outlined by Dr. Brau.

### Step 1

Place the patient in the supine position on an x-ray table.

The approach surgeon stands on the left and the assistant on the right. The level of the transverse incision in the craniocaudad plane depends on the level of the spine to be approached. A lateral x-ray of the spine is essential to determine the proper placement of this incision.




### **CAUTION**

If the patient is obese, avoid compressing the ulnar nerve when placing rail clamp.



#### **TIP**


When necessary, use a wider OR table or add  $2\frac{1}{4}$ " to the width of the table by using our Rail Extender (#41917).



### Step 2

The left rectus muscle is mobilized circumferentially. With the rectus muscle initially retracted medially, carefully incise the posterior sheath of transversalis fascia 4 to 5 cm until the peritoneum is seen to shine through. Grasp the edges with a hemostat and lift it away and very carefully dissect if from the peritoneum. Incise it as far inferiorly and superiorly as possible. Using your index finger, carefully push the peritoneum posteriorly at the edge of the fascial incision and slowly develop a plane between it and the undersurface of the internal oblique and transversus muscles and fascia. This will lead you to the retroperitoneal space.


Continue careful blunt finger dissection posteriorly, and then start pushing medially trying to elevate the peritoneum away from the psoas muscle. Be careful not to enter the retropsoas space at this point, as this will lead to unnecessary bleeding in a blind pouch. The genitofemoral nerve can be easily identified over the psoas. The ureter can usually be identified as the peritoneum is lifted away from the psoas. Both of these structures should be preserved from injury.



# Retroperitoneal Mini-Approach Set Up Steps (continued)


### Step 3

Once the psoas is identified, palpate medially to feel for the disc and vertebral body and iliac artery. At this point, if size of the incision allows, insert the entire hand and make a fist in the retroperitoneal area. Sweep with the closed fist up and down to elevate the peritoneum away in all directions. Continue with blunt dissection to expose the entire length of the common and external iliac arteries as far distally as possible, and then start careful blunt dissection along the lateral edge of the artery. This will expose the left common iliac vein just underneath the artery. Continue the dissection posteriorly to identify the ileolumbar vein(s). Variations in the formation of the common iliac vein and the lumbar veins are common, and great care must be exercised in order to identify, ligate and transect these veins and avoid avulsion. The left iliac vein and artery can now be separated away from the spine using gentle, peanut sponge, fingertip and blunt elevator dissection.




### Step 4

Secure the Elite II Rail Clamp to the table rail over the sterile drape on either side of the table, whichever side keeps the surgeon's operating field clear (A). The Crossbar is inserted into the joint and positioned 2 cm above the operative site (B).



### Step 5


All vascular structures are then swept from the left to right, providing adequate visualization of the disc(s) and vertebral bodies involved. Segmental vessels running across the valleys on the anterior surface of the bodies can be transected between clips and swept to the sides with blunt dissection. Make sure you can get at least one finger between the vein and the ligament so that you can palpate the right lateral edge of the spine with the vessels above your finger(s).



# Retroperitoneal Mini-Approach Set Up Steps (continued)


### Step 6

The lateral extension arms are attached to the crossbar (A) and positioned just above the horizon of the patient (B).




### Step 7

The surgeon's left hand then re-enters the retroperitoneal space with the rectus now moved laterally, and the fingers find their way to the right side of the spine. A Radiolucent Reverse Lip Anterior Spine Access blade is placed blindly on the right side of the spine using the finger(s) as a guide.



### Step 8


This blade is then attached to the lateral extension arms of the retractor frame, elevating the vascular structures and exposing the anterior surface of the spine.



Once secured to the Thompson Retractor, the reverse-lipped blade will not move. The reverse lip keeps the blade anchored to the edge of the spine and prevents it from slipping anteriorly once tension is applied. Without this reverse lip, the retractor blade will not work effectively.


### Step 9

With the rectus now retracted laterally, there will be much less resistance when pushing the retractor blade to expose the spine in a direct AP view and allow placement of the sleeves for insertion



of a threaded device, femoral ring or artificial disk. Place a second reverse-lip blade on the left side of the spine and attach to the Thompson frame. Commonly, additional retractor blades need to be placed superiorly and/or inferiorly to complete the exposure. With the blades well anchored to the lateral wall of the vertebral column, the spine surgeon and the assistant can now work on the disc without other hands or retractors being in the way and with relative security that vessels will not move around the retractors and expose themselves to injury.

# **Anterior Lumbar System Components**



### **ANTERIOR LUMBAR SYSTEM**

### #SL82019

REF	QTY	ITEM DESCRIPTION	PART #
А	1	Elite III Rail Clamp with 1 Cam Joint 16"	43902ACL
В	1	Crossbar with 2 Cam III Joints 24"	43900BC
С	2	20" Angled Arm (10" x 10") @ 45°	44120
D	2	24" Angled Arm (8" x 16") @ 45°	44124N
Е	5	Cam II Clip-on Quick Angle 8"	SL42126WGP
F	1	Micro-Adjustable II Clip-on Quick Angle 10"	SL45006CGP
G	2	T-Handle	60020
Н	1	Quick Angle Hand Held Adapter	SL42128G
I	1	Suction for Anterior Lumbar Surgery	51234
J	1	Anterior Lumbar Depth Gauge 335mm	51236
	2	Radio. Malleable Renal Vein 25mm x 140mm	SL46119BET
Κ	2	Radio. Malleable Renal Vein 25mm x 190mm	SL46119CET
L	2	Radio. Malleable 51mm x 203mm (2" x 8")	SL46122ET
М	1	Radio. Harrington 64mm x 152mm (2 ½" x 6")	SL46160ET
Ν	1	Radio. Splanchnich 51mm x 115mm (2" x 4 $\frac{1}{2}$ ")	SL46145ET
0	2	Radio. Balfour 65mm x 72mm (2 ½" x 2 ¾")	SL46140ET
	4	Radio. Concave 25mm x 100mm (1" x 4")	SL46560
	4	Radio. Concave 25mm x 120mm (1" x 4 ¾")	SL46562
Р	4	Radio. Concave 25mm x 140mm (1 " x 5 ½")	SL46564
	4	Radio. Concave 25mm x 160mm (1" x 6 1/4")	SL46566

REF	QTY	ITEM DESCRIPTION	PART #
	4	Radio. Concave 25mm x 180mm (1" x 7")	SL46568
Q	2	Reverse Lip Malleable 25mm x 200mm	SL46192ET
	2	Reverse Lip Malleable 32mm x 200mm	SL46193ET
	2	Reverse Lip Tapered 25mm x 100mm	SL46260TET
R	2	Reverse Lip Tapered 25mm x 150mm	SL46265TET
	2	Reverse Lip Tapered 25mm x 200mm	SL46270TET
	2	Reverse Lip Rigid 25mm x 110mm	SL46261RET
	2	Reverse Lip Rigid 25mm x 130mm	SL46263RET
S	2	Reverse Lip Rigid 25mm x 150mm	SL46265RET
	2	Reverse Lip Rigid 25mm x 170mm	SL46267RET
	2	Reverse Lip Rigid 25mm x 190mm	SL46269RET
	2	Reverse Lip Rigid 32mm x 120mm	SL46282RET
Т	2	Reverse Lip Rigid 32mm x 140mm	SL46284RET
	2	Reverse Lip Rigid 32mm x 160mm	SL46286RET
	2	Reverse Lip Rigid 32mm x 180mm	SL46288RET
	2	Reverse Lip Rigid 32mm x 200mm	SL46290RET
	1	Instrument Case	50000G
	1	Instrument Case - Exp. Reverse Lip ASA	50000ERL
	1	Instrument Case - AL Renal Vein Blades	50000ALR
	1	Instrument Case - Anterior Lumbar	50000ALT

# 0

### **SYSTEM NOTICE:**

Non-S-Lock system available, but may require a longer lead time. Please call for more information.

# **Table Adapters + Rail Extenders**

# Providing stable support on the OR table for the table mounted frame


When bedrail space is not available for a rail clamp, or, when the bedrail needs to be offset to accommodate obese patients, a rail adapter or rail extender should be applied.



### **TABLE ADAPTER**

Easily connects to Jackson Spine tables to add a standard bed rail for applying a table mounted Thompson Retractor System to your operation.

REF	ITEM DESCRIPTION	PART #
U	Jackson Spine Frame Adapter	41927



### **RAIL EXTENDERS**

Apply a Rail Extender to your OR table to increase your rail length or width and provide more attachment options for Thompson rail clamps.

REF	ITEM DESCRIPTION	PART #
V	Rail Extender 15" Single Clamp	5844
	Rail Extender 20" Single Clamp	41938
W	Rail Extender 14" Long with 2 ¼" Offset Dual Clamp	41917
	Rail Extender 22" Dual Clamp	41929

1.800.227.7543





phone: 231.922.0177 fax: 231.922.0174 thompsonsurgical.com

### EC REP Emergo Europe

Prinsessegracht 20 2514 AP The Hague THE NETHERLANDS

- \* Free trial valid for U.S. customers only. Customers outside U.S. please call +1-231-922-0177 for availability.
- © 2019 Thompson Surgical Instruments, Inc. Traverse City, Michigan. Printed in the U.S.A.
- ® S-Lock® Ø, PLA®, and the "T Circle" logomark 🏶 are Registered Trademarks of Thompson Surgical Instruments, Inc. Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064, US9872675 Other patents pending.

Symbol Legend:





