









Bariatric Liver Transplant



BARIATRIC SYSTEM #SL82001 Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43905AC	Elite III Rail Clamp w/2 Cam Joints 22"	2	Attaches to the rail of the OR table and holds the bilateral crossbar and angled arms in position.
B		41900PHL	41 1/2" Bilateral Crossbar Hinged (14" x 10" x 14")	1	Attaches to the rail clamps to provide bilateral stability, and holds the retractor handles in place around the incision.
C		44124W	24" Angled Arm (11" x 13") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
		50000BL	Instrument Case (26" x 10" x 5")	1	Holds instruments.
Tray 2					
D		SL42126WGP	Cam II Clip-on Quick Angle 8"	5	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down.
E		SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	2	Provides precise micro-adjustments.

BARIATRIC SYSTEM #SL82001 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		SL46119C	Concave Round Neck Renal 25mm x 191mm (1" x 7 1/2")	2	Renal vein or vessel retraction.
		SL46122A	Malleable 51mm x 254mm (2" x 10")	1	Soft tissue retraction.
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SO46124B	Malleable 102mm x 254mm (4" x 10")	1	
		SO46130	Kelly 76mm x 89mm (3" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SO46136A	Fence 102mm x 203mm (4" x 8")	1	Lateral retraction of abdominal wall on large procedures.
		SO46141	Balfour 83mm x 73mm (3 1/4" x 3")	2	Costal margin or abdominal wall retraction.
		SO46141D	Balfour 83mm x 123mm (3 1/4" x 5")	2	
		SO46141F	Balfour 83mm x 175mm (3 1/4" x 6 7/8")	1	
		SO46153	Richardson 51mm x 178mm (2" x 7")	1	Abdominal wall retraction.
		SL46161	Harrington 114mm x 254mm (4 1/2" x 10")	1	Organ or soft tissue retraction.
		SO46170C	Non-Slip Balfour R 77mm x 94mm (3" x 3 3/4")	1	Costal margin retraction. Aligns with anatomy, reducing slippage.
		SO46171C	Non-Slip Balfour L 77mm x 94mm (3" x 3 3/4")	1	
		SL46191	Malleable Finger 152mm (6")	1	Gentle organ retraction.
		SL46146	Splanchnic 51mm x 152mm (2" x 6")	1	Soft tissue retraction.
		50000G	Instrument Case (23"x 11" x 3 1/2")	1	Holds instruments.