



Abdominal Hysterectomy



Vaginal Hysterectomy



GYN SYSTEM #SL82010

Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43903AC	Elite III Rail Clamp with 2 Cam Joints 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	24" Crossbar with 2 Cam III Joints	1	Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles.
C		44124W	24" Angled Arm (11" x 13") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
D		44118N	18" Angled Arm (8" x 10") @ 45°	2	
E		44118C	Straight Arm 18" with 1 Cam Joint	1	Attaches to the rail clamp and holds hysterectomy tray in place.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
Tray 2					
F		41933	Hysterectomy Tray	1	Attaches to the straight arm and holds surgical instruments.
G		SL42126WGP	Cam II Clip-on Quick Angle 8"	3	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
H		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.

GYN SYSTEM #SL82010 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
J		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
K		SL46180	Heany 25mm x 102mm (1" x 4")	1	Abdominal or incisional wall retraction.
L		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
M		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	Soft tissue retraction.
N		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	
O		SL46118A	Malleable 25mm x 254mm (1" x 10")	2	
P		SL46118P	Malleable 25mm x 102mm (1" x 4")	2	Soft tissue retraction.
Q		SL46118	Malleable 25mm x 203mm (1" x 8")	2	
R		SL46122	Malleable 51mm x 203mm (2" x 8")	2	
S		SO46135	Fence 102mm x 127mm (4" x 5")	1	Lateral retraction of abdominal wall on large procedures.
T		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
U		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
V		5844	Rail Extender 15" with 1 Power Rail Clamp	1	Attaches to OR table rail to increase rail length or width and provide more attachment options for the rail clamps.