


HIPEC













HIPEC SYSTEM #SL82004

Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	2	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	24" Crossbar with 2 Cam III Joints (11" x 13")	1	Attaches to the rail clamp and holds the angled arms or lower abdominal bar in place around the incision. May also support retractor handles.
C		44014C	14" Extension Arm with Clip Joint	1	Attaches to the rail clamp and provides bilateral stability to the lower abdominal bar components in HIPEC procedures.
D		44122DB	22" Double Angled Arm	2	Attaches to the crossbar and angled closure arm to enclose the frame and holds the retractor handles in place around the incision.
		44112SP	20" Lower Abdominal Bar (14" x 6")	1	Attaches to the rail clamps or crossbar and holds the retractor handles in place around the incision.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.

**HIPEC SYSTEM #SL82004 (Continued)**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 2</b>					
E		44113B	13" Angled Closure Arm (6 1/2" x 6 1/2")	1	Attaches to the double angled arms, enclosing the frame around the incision.
		44118B	18" Angled Closure Arm (9" x 9")	1	
F		SL42126WGP	Cam II Clip-on Quick Angle 8"	6	Attaches retractor blades to the frame and angles retractor blades as needed, 45° up or down.
		SL42126WGLP	Cam II Clip-on Quick Angle 10"	2	
		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	6	Abdominal or incisional wall retraction.
		SO46140C	Balfour [sm lips] 66mm x 73mm (2 1/2" x 3")	2	Costal margin or abdominal wall retraction.
		SL46174	Deaver 51mm x 127mm (2" x 5")	2	Soft tissue retraction.
		SL46191	Malleable Finger 152mm (6")	1	Gentle organ retraction.
		SO46170B	Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2")	1	Costal margin retraction. Aligns with anatomy, reducing slippage.
		SO46171B	Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2")	1	
		SO46188	St Marks [with lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.