



Elite Rigid Holder + Flexible Holder



Elite Articulating Arm + Flexible Holder




LAPAROSCOPIC PRODUCTS FOR ELITE RAIL CLAMP


REF	IMAGE	PART #	ITEM DESCRIPTION	APPLICATION
A		43903AC	Elite III Rail Clamp with 2 Cam Joints 18"	Attaches to the rail of the OR table and holds the laparoscopic holding arm(s) in position.
B		42133DEN	Flexible Laparoscopic Holder for Elite	Unlimited range of motion and positioning with adjustable tension from flexible to rigid.
C		42139EU	Universal Laparoscopic Holder for Elite	Precise, rigid positioning for instruments requiring few adjustments or holding excess weight. Additional rigid holder options available.
		42139E	Rigid Laparoscopic Holder for Elite	
D		42138LE	Laparoscopic Articulating Arm for Elite	Holds firmly with a wide range of motion.
E		33181	Adjustable Scope Grip 5-6mm / 10-11mm	Holds scope firmly while protecting shaft. Fits scopes 5-6 or 10-11mm.
F		SL42161	S-Lock Nathanson Adapter (required for S-Lock Nathanson Hook use)	Allows quick, easy attachment of S-Lock Nathanson Hooks to laparoscopic holder universal end.
G		SL42201	S-Lock Nathanson Hook Medium 5mm (203mm x 167mm x 16mm)	Organ retraction.

Additional attachments available, not shown.








PRC Rigid Holder + Flexible Holder



PRC Articulating Arm + Flexible Holder



LAPAROSCOPIC PRODUCTS FOR POWER RAIL CLAMP

REF	IMAGE	PART #	ITEM DESCRIPTION	APPLICATION
A		41950	Power Rail Clamp (PRC)	Attaches to the rail of the OR table and holds the laparoscopic holding arm(s) in position.
B		42133DN	Flexible Laparoscopic Holder for PRC	Unlimited range of motion and positioning with adjustable tension from flexible to rigid.
C		42139U	Universal Laparoscopic Holder for PRC	Precise, rigid positioning for instruments requiring few adjustments or holding excess weight. Additional rigid holder options available.
		42139N	Rigid Laparoscopic Holder for PRC	
D		42138L	Laparoscopic Articulating Arm for PRC	Holds firmly with a wide range of motion.
E		33181	Adjustable Scope Grip 5-6 / 10-11mm	Holds scope firmly while protecting shaft. Fits scopes 5-6 or 10-11mm.
F		SL42161	S-Lock Nathanson Adapter (required for S-Lock Nathanson Hook use)	Allows quick, easy attachment of S-Lock Nathanson Hooks to laparoscopic holder universal end.
G		SL42201	S-Lock Nathanson Hook Medium 5mm (203mm x 167mm x 16mm)	Organ retraction.

Additional attachments available, not shown.