


Pediatric Abdominal



PEDIATRIC SYSTEM #SL82006 Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902IRC	Infinite Height Rail Clamp with 1 Cam III Joint	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	24" Crossbar with 2 Cam III Joints	1	Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles.
C		44120	20" Angled Arm (10" x 10")	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
		44014S	12" Angled Arm (6" x 6")	2	
D		SL42126WGP	Cam II Clip-on Quick Angle 8" *	5	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
		SL45011CGP	Mini Micro-Adjustable Clip-on Quick Angle 9"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case - 23" x 11" x 3 1/2"	1	Holds instruments.

# PEDIATRIC SYSTEM #SL82006 ( Continued )

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 2</b>					
		SL46115	Malleable 13mm x 64mm (½" x 2 ½")	2	Soft tissue retraction.
		SL46116	Malleable 13mm x 152mm (½" x 6")	2	
		SL46117	Malleable 13mm x 203mm (½" x 8")	2	
		SL46118P	Soft Malleable 25mm x 102mm (1" x 4")	2	
		SL46122P	Malleable 51mm x 102mm (2" x 4")	2	
		SL46191SP	Malleable Finger 76mm (3")	1	Gentle organ retraction.
		SL46191P	Malleable Finger 102mm (4")	1	
		SL46125	Richardson 19mm x 25mm (¾" x 1")	1	Muscle retraction for exposure with small or mini-incision.
		SL46125A	Richardson 19mm x 50mm (¾" x 2")	1	
		SL46126	Richardson 25mm x 32mm (1" x 1 ¼")	1	
		SL46127	Kelly 38mm x 51mm (1 ½" x 2")	1	Abdominal or incisional wall retraction.
		SL46128	Kelly 51mm x 64mm (2" x 2 ½")	1	
		SL46140P	Balfour 76mm x 51mm (3" x 2")	1	Costal margin or abdominal wall retraction.
		SL46142	Balfour 32mm x 32mm (1 ¼" x 1 ¼")	2	
		SL46142P	Balfour 51mm x 54mm (2" x 2")	1	
		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
		SL46162P	Harrington 38mm x 127mm (1 ½" x 5")	1	Organ or soft tissue retraction.
		SL46170	Deaver 25mm x 127mm (1" x 5")	1	Soft tissue retraction.
		50000G	Instrument Case (23" x 11" x 3 ½")	1	Holds instruments for sterilization, storage, and transportation.

\* For infants: Screw Down Handle 8" (SL42126NGP) recommended.