



Perineal Prostatectomy



Retropubic Prostatectomy



UROLOGY SYSTEM #SL82012

Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	24" Crossbar with 2 Cam III Joints	1	Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles.
C		44118N	18" Angled Arm (8" x 10") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
		44124N	24" Angled Arm (8" x 16") @ 45°	2	
D		SL42126WGP	Cam II Clip-on Quick Angle 8"	3	Attaches retractor blades to the crossbar or angled arm and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-On Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	Use on QA Handle to angle blade 45° up/down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
Tray 2					
		SL46116	Malleable 13mm x 152mm (1/2" x 6")	1	Soft tissue retraction.
		SL46122	Malleable 51mm x 203mm (2" x 8")	1	
		SL46123	Malleable 38mm x 203mm (1 1/2" x 8")	1	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	1	

UROLOGY SYSTEM #SL82012 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	2	Abdominal or incisional wall retraction.
		SO46136	Fence 102mm x 152mm (4" x 6")	1	Lateral retraction of abdominal wall on large procedures.
		SL46137	Holtgrewe 60mm x 152mm (2 3/8" x 6")	1	Malleable bladder retractor, designed to protect bladder during resection of prostate gland.
		SL46140P	Balfour [with lips] 76mm x 51mm (3" x 2")	2	Costal margin or abdominal wall retraction.
		SL46142P	Balfour [with lips] 51mm x 54mm (2" x 2")	2	
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SO46151	Richardson 51mm x 127mm (2" x 5")	1	Abdominal wall retraction.
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46163	Harrington 64mm x 203mm (2 1/2" x 8")	1	
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SL45753AM	Prostate Retractor 8mm (3/8")	1	Prostate retraction.
		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	Soft tissue retraction.
		SL46661	Double Angled 25mm x 64mm (1" x 2 1/2")	2	
		SL46662	Double Angled 13mm x 76mm (1/2" x 3")	2	
		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	
		SL46215	Murphy Blunt 25mm x 13mm (1" x 1/2")	1	Skin and shallow muscle retraction.
		50000G	Instrument Case (22" x 11" x 3 1/2")	1	Holds instruments.